



1 Create Encounter

Patient Encounter Creation: Enter Details

Patient Details


Patient ID*:
 Patient Name:
 Encounter/Visit ID:
 Date of Birth: February (20) | 12 | AMT
 Gender: Male Female Other Prefer not to say
 Type of Diabetes: Type 1 Type 2 Unknown
 Patient Dilator Status: Dilated (Hydraloc) Not Dilated (Non-Hydraloc)
Protocol Details
 Protocol Type: 2-Field
Camera Details
 Camera Model: Canon CR-2 AP



4 Track Available Results

Create Encounter

1	Patient Name/ID	Encounter/ID	Created	EyeArt/ID	Status
1	Smith, Ryan	123456	2021-Jun-09 14:39	92393	Completed
2	John Doe	123457	2021-Jun-09 14:39	92394	Completed
3	John Doe	123458	2021-Jun-09 14:39	92395	Completed
4	John Doe	123459	2021-Jun-09 14:39	92396	Completed
5	John Doe	123460	2021-Jun-09 14:39	92397	Completed
6	John Doe	123461	2021-Jun-09 14:39	92398	Completed
7	John Doe	123462	2021-Jun-09 14:39	92399	Completed
8	John Doe	123463	2021-Jun-09 14:39	92400	Completed
9	John Doe	123464	2021-Jun-09 14:39	92401	Completed
10	John Doe	123465	2021-Jun-09 14:39	92402	Completed
11	John Doe	123466	2021-Jun-09 14:39	92403	Completed
12	John Doe	123467	2021-Jun-09 14:39	92404	Completed
13	John Doe	123468	2021-Jun-09 14:39	92405	Completed
14	John Doe	123469	2021-Jun-09 14:39	92406	Completed
15	John Doe	123470	2021-Jun-09 14:39	92407	Completed
16	John Doe	123471	2021-Jun-09 14:39	92408	Completed



6 Create PDF Report

EYENUK **EyeArt**
Diabetic Retinopathy Analysis Report

Patient Information **General Information**

Patient ID: p001 Referring Location: Diab Clinic
 Patient Name: Ryan Sangoli Referring Provider: Dr. John Doe
 Date of Birth: 1977-Apr-05 EyeArt Control ID: 88178
 Gender: Male Dilator Status: Not Dilated
 Encounter ID: e001 Exam Analysis Date: 2020-Feb-25 14:29

EyeArt Diabetic Retinopathy (DR) Exam Result Summary

Negative for more than mild DR in both eyes. Retest in 12 months.

Right Eye Results **Left Eye Results**

mtmDR Result: Negative for more than mild DR mtmDR Result: Negative for more than mild DR
 vDR Result: Negative for vision-threatening DR vDR Result: Negative for vision-threatening DR

Refer to the EyeArt System User Manual for additional details or e-mail us at support@eyenuk.com.



EyeArt result “Ungradable”

- If encounter states “images of sufficient quality not present for all required fields,” retry fundus photography. Administer dilation drops before retrying, if necessary and possible.
- Check images for proper resolution, quality, positioning, and eye and field labels.



Error message: “Server is unreachable”

- Check your Internet connection and restart.
- If your computer requires a proxy to access the Internet, confirm the *EyeArt* client proxy server settings are correct. Please refer to the *EyeArt* System User Manual for more information.



Error message: “Incorrect login credentials”

- Check your username and password.
- Check the license key, make sure there are no beginning or trailing spaces.
- Check your Internet connection.
- If you are using the correct credentials and still unable to login/register, please contact support@eyenuk.com.



Images do not appear in the *EyeArt* client capture window

- Check if the capture folder is set correctly in the *EyeArt* client capture window. If not, use the “Reset Capture Folder” button to change the capture folder.
- Confirm the camera is correctly configured to save images to the capture folder.

Refer to the *EyeArt* System User Manual for additional details or e-mail us at support@eyenuk.com.